


TWO MWA PERCUTANEOUS ANTENNAS

The following illustration and table depict the approximate size of an ablation zone that may be created by two Evident™ MW Ablation Percutaneous Antennas energized simultaneously at various power and time settings using various spacings.



S	30 W 1.5 cm	10:00 2.0 cm	45 W 1.5 cm	5:00 2.0 cm	45 W 1.5 cm	10:00 2.0 cm
A	3.9 cm	4.0 cm	3.6 cm	4.2 cm	4.2 cm	4.5 cm
B	4.1 cm	4.1 cm	4.1 cm	4.0 cm	4.5 cm	4.5 cm
C	3.0 cm	2.7 cm	2.8 cm	2.5 cm	3.5 cm	3.5 cm
M	0.6 cm	0.6 cm	0.6 cm	0.6 cm	0.7 cm	0.7 cm
	25 cm ³	23 cm ³	21 cm ³	22 cm ³	35 cm ³	37 cm ³