


# ELECTROSURGICAL ACCESSORIES

## EDGE™ Coated Electrodes


*The Clear Cut Solution to A Sticky Problem*

Safety Sleeve™ insulated electrodes provide an increased level of safety, especially in procedures where the pencil tip can come into contact with unintended tissue.

There is no need to insulate these electrodes with rubber catheters or electrode tip protectors.



Surgeons find EDGE™ coated blades to be equal or superior to other coated or "nonstick" blades across a wide variety of procedures.



EDGE™ coated blades have been proven to provide performance beyond all other coated blade electrodes on the market. EDGE™ needle electrodes are designed to deliver uncompromised cutting performance for fine dissection, no matter what type of tissue.

### EASY TO CLEAN

With Valleylab's proprietary coating, no time-consuming scraping is required. Just wipe the electrodes clean with a wet sponge.

### DURABLE AND BENDABLE

EDGE™ electrodes are designed to be flexible without compromising the integrity of the coating.

### UNCOMPROMISED CUT PERFORMANCE

No matter what type of tissue or procedure, EDGE™ electrodes provide cut performance like stainless steel – but with a clean advantage.

### CHOICE OF STYLES

Purchase EDGE™ electrodes alone, in your favorite Valleylab pencil, or as part of the AccuVac™ smoke evacuation pencil.

# EDGE™ Coated Electrodes

## Order Information

Specify  
**EDGE™ Coated Electrodes**  
and pencils in your  
surgical kit packs.

CATALOG NUMBER	DESCRIPTION	ORDER QUANTITY	CATALOG NUMBER	DESCRIPTION	ORDER QUANTITY
Sterile, single use electrodes			E1465	7.21 cm (2.84 in.) PTFE insulated coated needle electrode for use with all Valleylab pencils	50/case
E1450X	6.35 cm (2.5 in.) coated blade electrode for use with all Valleylab hex-locking pencils	50/case	E1465B	7.21 cm (2.84 in.) PTFE insulated coated needle electrode with Safety Sleeve™ for use with all Valleylab handswitching pencils*	25/case
E1450G	7.62 cm (3.0 in.) coated blade electrode for use with all nonhex-locking pencils	50/case	E1465-4	10.16 cm (4 in.) extended PTFE insulated coated needle electrode for use with all Valleylab pencils	50/case
E1450-4	10.16 cm (4 in.) extended coated blade electrode for use with all Valleylab pencils	50/case	E1465B-4	10.16 cm (4 in.) extended PTFE insulated coated needle electrode with Safety Sleeve™ for use with all Valleylab handswitching pencils*	25/case
E1450-6	16.51 cm (6.5 in.) extended coated blade electrode for use with all Valleylab pencils	50/case	E1465-6	16.51 cm (6.5 in.) extended PTFE insulated coated needle electrode for use with all Valleylab pencils	50/case
E1452	7.21 cm (2.84 in.) coated needle electrode for use with all Valleylab pencils	50/case	E1475X	6.99 cm (2.75 in.) coated blade electrode for use with all Valleylab hex-locking pencils.	50/case
E1452-6	16.51 cm (6.5 in.) extended coated needle electrode for use with all Valleylab pencils	50/case	Sterile, single use electrosurgical pencils with EDGE™ Coated Blade		
E1455	6.99 cm (2.75 in.) PTFE insulated coated blade electrode for use with all Valleylab pencils	50/case	E2250H	Footswitching pencil with holster and 3 m (10 ft) cord	50/case
E1455B	6.99 cm (2.75 in.) PTFE insulated coated blade electrode with Safety Sleeve™ for use with all Valleylab handswitching pencils*	25/case	E2350H	Rockerswitch pencil with holster and 3 m (10 ft) cord	50/case
E1455-4	10.16 cm (4 in.) extended PTFE insulated coated blade electrode for use with all Valleylab pencils	50/case	E2350HS	Rockerswitch pencil with AccuVac™ smoke evacuation attachment	25/case
E1455B-4	10.16 cm (4 in.) extended PTFE insulated coated blade electrode with Safety Sleeve™ for use with all Valleylab handswitching pencils*	25/case	E2450H	Buttonswitch pencil with holster and 3 m (10 ft) cord	50/case
E1455-6	16.51 cm (6.5 in.) extended PTFE insulated coated blade electrode for use with all Valleylab pencils	50/case	E2450HS	Buttonswitch pencil with AccuVac™ smoke evacuation attachment	25/case

E2250H

E2350H

E2450H



**Valleylab**  
a division of Tyco Healthcare Group LP  
5920 Longbow Drive  
Boulder, CO 80301-3299 USA  
Valleylab Direct: 800 255 8522  
To Order Product: 800 722 8772

**tyco** / Healthcare