

Valleylab Footswitch Controls

Monopolar and bipolar power at your command



STANDARD BIPOLAR FOOTSWITCH (E6009)

allows for independent control of bipolar output with all Valleylab Force series generators.

DOME STYLE BIPOLAR FOOTSWITCH (E6019)

may be activated from any angle and has a 15 foot cord. Allows for independent control of bipolar output with all Valleylab Force series generators.

NEW LOW PROFILE MONOPOLAR FOOTSWITCH (E6008B)

is designed specifically for Valleylab full powered generators. Manufactured of high quality cast aluminum and features highly visible color coded foot pedals for safety. It also features a non-skid rubber base, water resistant strain relief connectors, and a divider bar between the CUT and COAG pedals.

MONOPOLAR FOOTSWITCH (E6008)


is designed for Valleylab full powered generators. Safety features include a water resistant design, a divider bar between the CUT and COAG pedals, and a non-skid rubber base. It also includes a water-resistant strain relief connector.

HIGH QUALITY

makes Valleylab the world leader in electrosurgical technology. These durable footswitches are manufactured to the highest standards of quality for exceptional performance and reliability during surgery.

Footswitch Controls

Technical Specifications



E6006-11

MONOPOLAR FOOTSWITCH

designed for use with the Valleylab SurgiStat™ electro-surgical generator only.

WEIGHT AND DIMENSIONS

	E6008B Low Profile Monopolar Footswitch	E6008 Standard Monopolar Footswitch	E6009 Standard Bipolar Footswitch	E6019 Dome Bipolar Footswitch
Length	19.1 cm (7.5 in.)	25.4 cm (10 in.)	11.7 cm (4.6 in.)	n/a
Width	21.8 cm (8.6 in.)	30.5 cm (12 in.)	8.9 cm (3.5 in.)	9.5 cm (3.75 in.)
Height	5.7 cm (2.25 in.)	14.2 cm (5.6 in.)	4.6 cm (1.8 in.)	1.3 in.
Weight	1.86 kg (4.1 lbs.)	2.95 kg (6.5 lbs.)	0.77 kg (1.7 lbs.)	0.77 kg (1.7 lbs.)
Cord Length	3.0 m (10 ft.)	3.0 m (10 ft.)	4.6 m (15 ft.)	4.6 m (15 ft.)

ORDER INFORMATION

CATALOG NUMBER	DESCRIPTION	ORDER QUANTITY
E6008B	Low profile compact monopolar footswitch for use with all full powered Valleylab electro-surgical generators.	1 each
E6008	Standard monopolar footswitch for use with all full powered Valleylab electro-surgical generators.	1 each
E6009	Bipolar footswitch designed for independent control of bipolar output with all Valleylab Force series electro-surgical generators.	1 each
E6019	Dome style bipolar footswitch for independent control of bipolar Valleylab Force series generators.	1 each
E6006-11	Monopolar footswitch designed for the Valleylab low powered SurgiStat™ electro-surgical generator only.	1 each

World Headquarters
Valleylab
5920 Longbow Drive
Boulder, CO 80301-3299 USA
800 255 8522



An ISO 9001 certified company