

Opti²™ Electrodes

Optimum Performance for Laparoscopic Procedures



CONVENIENT

use with your Valleylab handswitching pencils ensures the maximum in efficiency as well as cut and coag capabilities.

FLEXIBLE CONFIGURATION

options allow you to select the electrode and handswitching pencil that best suit your specific surgical need. Opti²™ single use solid electrodes are also interchangeable in the sterile field.

SAFETY AND PERFORMANCE

are integral to the design of the insulated, 28 cm shaft, solid electrodes. The Opti²™ single use electrodes feature a 4.8 mm diameter to fit in 5 mm ports and Valleylab's unique anti-rotation hex-lock fitting which secures the electrode in the handset throughout the procedure.

Handset and Opti²™ Single Use Solid Electrodes

Technical Specifications

OPTI²™ SINGLE USE SOLID LAPAROSCOPIC ELECTRODES

| | |
|-----------------|--|
| Working Length | 11 in. (28 cm) |
| Diameter | 4.8 mm |
| Insulation | Polyolefin |
| Compatible with | Valleylab E2515, E2516, and E2517H pencils |

HANDSWITCHING PENCILS

| | |
|--------------------------------------|--|
| Electrosurgical Cut and Coag Control | Rocker switch (E2515), Button switch (E2516), Button switch with power control (E2517H) |
| Cord | Three conductor, 10 ft. (3.0 m) length, terminated with plug for direct use with generator |
| Compatible with | All Valleylab generators |

ORDER INFORMATION

| CATALOG NUMBER | DESCRIPTION | ORDER QUANTITY |
|--|------------------------------|----------------|
| Sterile Solid Laparoscopic Electrodes | | |
| E2770-28 | Straight spatula tip | 10/case |
| E2771-28 | Curved spatula tip | 10/case |
| E2772-28 | Wire J hook tip | 10/case |
| E2773-28 | Wire L hook tip | 10/case |
| E2774-28 | Flat L hook tip | 10/case |
| E2775-28 | Ball tip | 10/case |
| Sterile Electrosurgical Pencils | | |
| E2515 | Handswitching, rocker switch | 50/case |
| E2516 | Handswitching, button switch | 50/case |
| E2517H | Handswitching, power control | 10/case |



E2770-28



E2771-28



E2772-28



E2773-28



E2774-28



E2775-28

World Headquarters
Valleylab Inc
a division of United States Surgical Corporation
5920 Longbow Drive
Boulder, CO 80301-3299 USA
800 255 8522



An ISO 9001 certified company