

ELECTROSURGICAL LAPAROSCOPIC INSTRUMENTATION

Opti⁴_{TM} Electrodes and Handset

Optimum Performance for Laparoscopic Procedures

E2781R-28ASP

E2784R-28ASP

E2783R-28ASP

E2782R-28ASP

E2750

ENHANCED FEATURES

Aspiration holes and additional electrode tip insulation aid in laparoscopic procedures requiring precise dissection with decreased risk of collateral spread and increased suction demands.

FLEXIBLE DESIGN

Options include a choice between retractable and nonretractable electrodes, pencil grip handset, and multiple tip designs. Optimize specific surgical demands safely and easily with the Opti⁴_{TM} electrodes.

FULL PERFORMANCE

Capabilities include suction, irrigation, cutting and coagulation, and electrode retractability incorporated into one instrument. The Opti⁴_{TM} multi-function device increases efficiency by reducing the number of instruments required.

Laparoscopic Handset and Electrodes

Technical Specifications



E2781R-28ASP*



E2782R-28ASP*



E2783R-28ASP*



E2784R-28ASP*



E2780-28



E2786-28

OPTI⁴™ SINGLE USE LAPAROSCOPIC HANDSET

Configuration	Pencil
Suction/Irrigation Control	Button activated
Electrosurgical Cut and Coag Control	Rocker switch activated
Suction Tube Fitting	4.7 mm (3/16 in.) barb connector
Irrigation Tube Fitting	Female locking luer connector
Cord	Attached three conductor, 3.0 m (10 ft.) length, terminated with plug for direct use with generator
Compatible with	All Valleylab generators and E2781R-E2784R series retractable laparoscopic electrodes, and E2780 series nonretractable laparoscopic electrodes.

OPTI⁴™ SINGLE USE LAPAROSCOPIC ELECTRODES

Working Length	28 cm (11 in.) and 36 cm (14 in.)
Diameter	4.8 mm
Insulation	Polyolefin (E2780 series) Polycarbonate and PTFE (E2781R-E2784R series)
Compatible with	Valleylab E2750 laparoscopic handset

Order Information

CATALOG NUMBER	DESCRIPTION	ORDER QUANTITY
Sterile Laparoscopic Handset		
E2750	Pencil grip handset	10/case
Sterile, Nonretractable Laparoscopic Electrodes		
E2780-28	Straight spatula tip, 28 cm	10/case
E2781-28	Curved spatula tip, 28 cm	10/case
E2782-28	Wire J hook tip, 28 cm	10/case
E2783-28	Wire L hook tip, 28 cm	10/case
E2784-28	Flat L hook tip, 28 cm	10/case
E2786-28	Cylindrical tip, 28 cm	10/case
Sterile, Retractable Laparoscopic Electrodes, with aspiration holes		
E2781R-28ASP	Curved spatula tip, 28 cm	10/case
E2782R-28ASP	Wire J hook tip, 28 cm	10/case
E2783R-28ASP	Wire L hook tip, 28 cm	10/case
E2783R-36ASP	Wire L hook tip, 36 cm	10/case
E2784R-28ASP	Flat L hook tip, 28 cm	10/case

Valleylab
a division of Tyco Healthcare Group LP
5920 Longbow Drive
Boulder, CO 80301-3299 USA
800 255 8522

© 2004 Valleylab

Specifications subject to change without notice.

tyco / Healthcare / **Valleylab**

Printed in USA

945605402

11/04